

STERILE EXAMINATION GLOVES

Class I Sterile / Cat. I PPE



proHAND®

Material	Latex	Size Ratio	S - XL
Force At Break /N/	BA Min. 6 / AA Min. 6	Color	Milky
AQL	1.0	Surface	Smooth Powdered
Thickness	F: 0,11 P: 0,10 C: 0,08 mm	Pack Type	4 x 50
Lenght	Min. 240 mm	Application	First Aid Rooms, Dental Surgery, Diagnostics, Nursing Care



proHAND® PF NITRILE

Material	Nitrile	Size Ratio	S - XL
Force At Break /N/	BA Min. 6 / AA Min. 6	Color	Blue
AQL	1.0	Surface	Finger Textured Chlorinated
Thickness	F: 0,05 P: 0,05 C: n/a mm	Pack type	4 x 50
Lenght	Min. 240 mm	Application	First Aid Rooms, Dental Surgery, Diagnostics, Nursing Care