

GAZE COMPRESSE STERILE



MAIN SPECIFICATIONS

- Totally free of loose threads, even when planned the first fold.
- Sterilized.
- Presentation: Disposable bags type "peel pack".
- Sterilization batch, expiration date and sterilization indicator in each package.
- Composition: 100% cotton wool, according to European Standard 14079.
- Available in various dimensions, folds, layers, wire numbers / cm², wired or wireless contrast X-ray and with or without security tape.
- Packaging mixed PA / PP and paper.
- Packaging single or double for added security
- The compresses are produced in order to meet the European regulatory requirements for classified according to ISO standards environment rooms.

Code	Description	Selling unit
01043010	GAZE COMPRESSE STERILE 5X5 CM 8 LAYERS 17 WIRES/CM ²	CX.100
01043011	GAZE COMPRESSE STERILE 7,5X7,5 CM 8 LAYERS 17 WIRES/CM ²	CX.100

01043060	GAZE COMPRESSE STERILE 10X10 CM 10 LAYERS 17 WIRES/CM2	CX.100
01043017	GAZE COMPRESSE STERILE 10X20 CM 8 LAYERS 17 WIRES/CM2	CX.360
01043056	GAZE COMPRESSE STERILE 7,5X7,5 CM	CX.100
01043000	GAZE COMPRESSE STERILE 10X10 CM PACKAGE 5	CX. 200
01043018	GAZE COMPRESSE STERILE 40X50 CM PACKAGE 3	CX. 108
01043019	GAZE COMPRESSE STERILE 20X15 CM PACKAGE 5	CX. 350
01043014	GAZE COMPRESSE STERILE CONTRAST X-RAY 20X15 CM PACKAGE 5	CX. 80
01043004	GAZE COMPRESSE STERILE CONTRAST X-RAY 40X50 CM PACKAGE 5	CX.110
01043038	GAZE COMPRESSE STERILE CONTRAST X-RAY 10X10 CM PACKAGE 5	CX. 200
01043053	GAZE COMPRESSE STERILE CONTRAST X-RAY 10X10 CM PACKAGE 10	CX. 600
01043021	GAZE COMPRESSE STERILE 10X10 CM PACKAGE 10	CX.100
01043022	GAZE COMPRESSE STERILE 10X10 CM PACKAGE 5	CX.200
01043023	GAZE COMPRESSE STERILE 10X20 CM PACKAGE 5	CX.100
01043025	GAZE COMPRESSE STERILE 5X5 CM PACKAGE 5	CX. 150
01043027	GAZE COMPRESSE STERILE 10X10 CM PACKAGE 10	CX. 20
01043036	GAZE COMPRESSE STERILE 7,5X7,5 CM PACKAGE 5	CX. 330
01043039	GAZE COMPRESSE STERILE 7,5X7,5 CM PACKAGE 5	CX.280
01043037	GAZE COMPRESSE STERILE 20X40 CM PACKAGE 5	CX. 30
01043040	GAZE COMPRESSE STERILE 20X40 CM PACKAGE 5	CX. 20
01043041	GAZE COMPRESSE STERILE 5X5 CM PACKAGE 5	CX. 36
01043031	GAZE COMPRESSE STERILE 20X20 CM PACKAGE 5	CX. 60
01043052	GAZE COMPRESSE STERILE CONTRAST X-RAY 10X10 CM 8 LAYERS 17 WIRES/CM2	CX.200
01043042	GAZE COMPRESSE STERILE 10X10 CM 8 LAYERS 17 WIRES/CM2	CX.100