

## ALTEADORES FOR TOILET SEATS



### MAIN SPECIFICATIONS

- Made of high density polyethylene.
- Solid and sturdy structure.
- Designed anatomically to give a better pelvic positioning and allow good distribution of body weight over the entire contact area.
- Extremely safe with adjustable closure to any shape and size toilet.
- They contain anti-slip rubber on each side, ensuring greater lateral stability.
- Front concavity to facilitate personal hygiene.
- With handling instructions, assembly and cleaning included in individual packaging.
- Available in two models with different heights: Simple model (without arms), 14 cm and Model arms with 10 cm high.
- The model with arms is extremely useful for the user can be up without any help from others.
- Arms aluminum lined anti-slip sponge for added security.

Code	Description	Selling unit
31009002	ALTEADOR TOILET SEAT WITH ADJUSTABLE ARMS SUPPORT	CX. 1