

HYDROFILM ROLL TRANSPARENT FILM



MAIN SPECIFICATIONS

- Hydrofilm roll is a film dressing for transparent fixation with varied uses. It offers protection against dirt and germs and can be cut exactly to size. For all-over dressing cover and fixation, especially when you want to monitor exudate absorption. Also as a waterproof dressing cover and dressing retention on joints, round and conical body parts.
- Waterproof film dressing made of transparent, breathable polyurethane film, coated with hypoallergenic adhesive; can be unrolled and easily cut to size thanks to its backing paper which aids in the application; no maceration, adheres securely, can be removed painlessly and without leaving any residues; sterilizable by ethylene oxide or radiation, indifferent to radiation.

Code	Description	Selling unit
03046027	HYDROFILM ROLL 10 M X 5 CM	CX. 1
03046028	HYDROFILM ROLL 10 M X 10 CM	CX. 1