

# HYDROFILM ROLL TRANSPARENT FILM



## MAIN SPECIFICATIONS

- Hydrofilm roll is a film dressing for transparent fixation with varied uses. It offers protection against dirt and germs and can be cut exactly to size. For all-over dressing cover and fixation, especially when you want to monitor exudate absorption. Also as a waterproof dressing cover and dressing retention on joints, round and conical body parts.
- Waterproof film dressing made of transparent, breathable polyurethane film, coated with hypoallergenic adhesive; can be unrolled and easily cut to size thanks to its backing paper which aids in the application; no maceration, adheres securely, can be removed painlessly and without leaving any residues; sterilizable by ethylene oxide or radiation, indifferent to radiation.

| Code     | Description                 | Selling unit |
|----------|-----------------------------|--------------|
| 03046027 | HYDROFILM ROLL 10 M X 5 CM  | CX. 1        |
| 03046028 | HYDROFILM ROLL 10 M X 10 CM | CX. 1        |