

LEUKOSILK SILK ADHESIVE



MAIN SPECIFICATIONS

- The ideal tape for sensitive skin
- Leukosilk is made of a fine acetate fabric coated with a skin-friendly polyacrylate adhesive.
- The product conforms well to body contours and is gentle on the skin. It can be torn easily by hand without fraying of the edges. Leukosilk is universal in use.
- Main uses: Fixation of wound dressings of all types and sizes on sensitive skin, securing of tubes, catheters, probes and cannulae.

Code	Description	Selling unit
3254024	LEUKOSILK ADHESIVE 5 M X 2,5 CM	CX. 12
3254025	LEUKOSILK ADHESIVE 5 M X 5 CM	CX. 6
3254023	LEUKOSILK ADHESIVE 5 M X 1,25 CM	CX. 24