

Negatoscope



MAIN SPECIFICATIONS

- Steel sheet structure;
- Epoxy paint;
- Front panel in white opal acrylic glass;
- Flushing electricity with 1 circular lamp;
- With earthing plug;
- Stainless steel plate for x-ray fixation / radiographs.
- Weight: 6 Kg
- Dimensions: 115 x 375 x 450mm

Code	Description	Selling unit
20019001	NEGATOSCÓPIO SIMPLES (1900)	UN

FACTORMED
Produtos Médicos & Hospitalares