

## SURGICAL GREEN MASKS



### MAIN SPECIFICATIONS

- Made of non-woven material in green.
- Latex-free, fibreglass-free.
- Hypoallergenic.
- Mask properties 3-ply mask with 2 round earloops or tie-on.
- Shingle pleat design prevents fluid pooling and cross-contamination.
- Adjustable nosepiece to customise your own desired fit.
- Single use.

Code	Description	Selling unit
13011040	SURGICAL GREEN MASK TRIPLE LAYERS WITH TIE ON	CX.50
13008008	SURGICAL GREEN MASK TRIPLE LAYERS WITH ELASTICS	CX.50