

Marquesa ginecológica e obstetrícia



MAIN SPECIFICATIONS

- Structure in rectangular steel tube;
- With bed in MDF and nappa-coated foam;
- Divided into 3 sections, 2 detachable and 1 articulated;
- Manual overhead section for gas shock absorber system;
- Adjustable height for electric motor with foot pedal and hand control for trendelenburg position;
- Epoxy paint;
- Base with 4 wheels, 2 with brake;
- Size: 1980x660x (600 to 1020) mm
- Weight: 80 kg

Code	Description	Selling unit
24009003	MARQUESA OBSERVAÇÃO E TRATAMENTO GINECOLÓGICO E OBSTETRÍCIO	UN