

GRAB RAIL FOR SETTING TO THE WALL FOR WC



MAIN SPECIFICATIONS

- With 2 Structure in round stainless steel tube 35 mm;
- With 3 points of attachment to the wall;
- With 2 angles of 90°;
- Dimensions 610x300 mm
- Weight: 2.3 kg

Code	Description	Selling unit
AP 4116	GRAB RAIL FOR SETTING TO THE WALL FOR WC LEFT	UN
AP 4115	GRAB RAIL FOR SETTING TO THE WALL FOR WC RIGHT	UN