

OMNIPLAST ADHESIVE



MAIN SPECIFICATIONS

- Hypoallergenic adhesive textile fabric adhesive synthetic glue.
- Resistant to water, but permeable to air and water vapor, adhering securely.
- Leaves no residue on the skin.
- Easy to tear due to its banks in zig-zag.
- It does not absorb X-rays; weather resistant and temperature resistant.
- For all the variety of dressings retention to fix probes, cannulas, etc.
- For patients with normal skin.

| Code | Description | Selling unit |
|----------|-----------------------------------|--------------|
| 03046006 | OMNIPLAST ADHESIVE 5 M X1,25 CM | CX. 20 |
| 03046007 | OMNIPLAST ADHESIVE 5 M X 2,5 CM | CX. 20 |
| 03046008 | OMNIPLAST ADHESIVE 5 M X 5 CM | CX. 20 |
| 03046009 | OMNIPLAST ADHESIVE 9,2 M X1,25 CM | CX. 24 |
| 03046010 | OMNIPLAST ADHESIVE 9,2 M X 2,5 CM | CX.12 |

| | | |
|----------|--------------------------------|------|
| 03046011 | ADESIVO OMNIPLAST 9,2 M X 5 CM | CX.6 |
|----------|--------------------------------|------|

| | | |
|----------|--------------------------------------|------|
| 03046015 | OMNIPLAST ADHESIVE WHITE 10 M X10 CM | CX.5 |
|----------|--------------------------------------|------|
