

OMNIPLAST ADHESIVE



MAIN SPECIFICATIONS

- Hypoallergenic adhesive textile fabric adhesive synthetic glue.
- Resistant to water, but permeable to air and water vapor, adhering securely.
- Leaves no residue on the skin.
- Easy to tear due to its banks in zig-zag.
- It does not absorb X-rays; weather resistant and temperature resistant.
- For all the variety of dressings retention to fix probes, cannulas, etc.
- For patients with normal skin.

Code	Description	Selling unit
03046006	OMNIPLAST ADHESIVE 5 M X1,25 CM	CX. 20
03046007	OMNIPLAST ADHESIVE 5 M X 2,5 CM	CX. 20
03046008	OMNIPLAST ADHESIVE 5 M X 5 CM	CX. 20
03046009	OMNIPLAST ADHESIVE 9,2 M X1,25 CM	CX. 24
03046010	OMNIPLAST ADHESIVE 9,2 M X 2,5 CM	CX.12

03046011	ADESIVO OMNIPLAST 9,2 M X 5 CM	CX.6
----------	--------------------------------	------

03046015	OMNIPLAST ADHESIVE WHITE 10 M X10 CM	CX.5
----------	--------------------------------------	------
