

EYCOPAD EYEPIECE DRESSING STERILE



MAIN SPECIFICATIONS

- Ocular dressing with excellent cushion effect.
- Dressing eye made of cotton, gauze layer of cotton wool.
- Particularly soft and absorbent.
- Excellent cushion effect which is maintained for long periods of time.
- Sterilized and packed individually.
- For the treatment of traumatic injuries of the eye and post treatment operations to the eye.

Code	Description	Selling unit
01046008	EYEPIECE STERILE DRESSING EYCOPAD 56X70 MM	CX.25