

Folding triple screen panel



MAIN SPECIFICATIONS

- Structure in round steel tube;
- Epoxy and Stainless Steel Paint
- 3 panels in white plastic screen;
- Size: 680x1685 mm (measures of panel)

Code	Description	Selling unit
24019022	FOLDING TRIPLE SCREEN PANEL EPOXY WITH WHEELS	UN
BI 1211	FOLDING TRIPLE SCREEN PANEL EPOXY WITHOUT WHEELS	UN
BI1212	FOLDING TRIPLE SCREEN PANEL STAINLESS STEEL WITH WHEELS	UN
BI 1213	FOLDING TRIPLE SCREEN PANEL INOX WITHOUT WHEELS	