

ROTATING BATH SEAT



MAIN SPECIFICATIONS

- Structure in stainless steel tube;
- Fixed to the floor (by means of bushes and screws);
- PVC seat and backrest;
- 2 wrists at the ends of the arms;
- Elevation by means of spindle (crank);
- Rotary movement at 180°.

| Code | Description | Selling unit |
|----------|--|--------------|
| CD 1305 | ROTATING BATH SEAT | UN |
| 31009004 | CADEIRA DE BANHO FIXA COM RODAS 44 CM CAD6 | UN |