

## MOBILE SHOWER CHAIR



### MAIN SPECIFICATIONS

- Structure in stainless steel tube;
- With 4 wheels, 2 with brake;
- PVC seat and backrest;
- 4 wrists at the ends of the arms.
- Adjustable in 3 heights
- Dimensions: 600 x 600 x 900 | (450 to 580) mm
- Weight: 18 kg

Code	Description	Selling unit
31009022	CADEIRA DE BANHO C/ 2 RODIZIOS + RODA GRANDE MACIÇA CAD3	UN
31009023	CADEIRA DE BANHO C/ 2 RODIZIOS + RODA GRANDE MACIÇA CAD5	UN
31009018	CADEIRA DE BANHO C/ 4 RODIZIOS CAD3	UN
CD 1300	MOBILE SHOWER CHAIR	UN