

## SHOWER TROLLEY



### MAIN SPECIFICATIONS

- Steel tube structure, epoxy paint
- Bed in phenolic compact, tub in vulcanized fabric, with water flow through valve
- 4 folding rails, electrical or hydraulic lift by pedal, train position and anti-train driven by 2 gas springs
- 4 wheels being 2 with brake
- Dimensions: 1860 x 700 x (430 to 850) mm
- Weight: 73 kg
- Height adjustable
- Overall Overall Width 7600 mm

Code	Description	Selling unit
24019041	MACA DE BANHO HIDRÁULICA	UN
24019042	MACA DE BANHO ELÉCTRICA	UN