

DISPOSABLE MANUAL MASK



MAIN SPECIFICATIONS

- Mask Pre crowded - The block is inflated to an optimum level to ensure maximum grip and reduce the risk of leakage. Eliminates the risk of allergic reactions to latex, the mask does not contain latex.
- Eliminates the risk of cross-contamination - single-use masks, do not produce a single cross-contamination.
- Transparent like glass, allows visually check the condition of the patient and provides good adhesion.
- ISO Connectors - Ensures compatibility with all connectors and breathing systems with standard connectors of 22 mm and 15 mm incorporating.
- Sizes - available in various sizes, from newborn to large adult size.

Code	Description	Selling unit
39011040	DISPOSABLE MANUAL MASK NR 0	CX. 1
39011041	DISPOSABLE MANUAL MASK NR 1	CX. 1
39011042	DISPOSABLE MANUAL MASK NR 2	CX. 1
39011043	DISPOSABLE MANUAL MASK NR 3	CX. 1

39011044	DISPOSABLE MANUAL MASK NR 4	CX. 1
39011045	DISPOSABLE MANUAL MASK NR 5	CX. 1
39011038	DISPOSABLE MANUAL MASK NR 00	CX. 1