

# OVERBOOTS



## MAIN SPECIFICATIONS

- Composition: Non woven PLP SPUN 40gr/sqm + white PE sole 80?
- Color White PLP & White sole
- Thickness/Weight: 32 grs/pair
- Size: 370x325 mm
- Efficient anti-sliding power
- The height of product covers the lower limbs to the bottom of the calf. This product can be used with a protective overall
- Thanks to high resistance of sole, limited linting
- Good resistance to blood, ethilic alcohol, ether, acetone, ink, resistant to ethylen Oxide sterilization and radiation

Code	Description	Selling unit
13008021	TRANSPARENT BOOT COVERS 500 PIECES	CX. 500