

## ADUL BIBS WITH POCKET



### MAIN SPECIFICATIONS

- For patient and clothes protection from fluids and splashes
- Material Made of 1-ply absorbent tissue (20 g) and 1-ply PE (13 g) for ultimate protection.
- It has an embossed pattern.
- Composition of tissue: 100% pulp
- Type Adult bibs are designed with a flip over pocket that maintains its shape during usage.
- Soft and strong upper bib part with reinforced 3-ply non-woven strip.
- Non-toxic, non-allergenic, non-sterile
- Size Bib size: 66 x 37 cm
- Pocket size: 10 x 37 cm

Code	Description	Selling unit
13008001	ADULT BIBS WITH POCKET 37X67CM	CX.100