

RESUSCITATOR WITH MASK



MAIN SPECIFICATIONS

- The anatomical design of the masks allows for better comfort in the face of the patient.
- The valve can rotate 360 degrees, allowing the user to use both hands.
- Not autoclavable, synthetic, disinfected by liquids.
- All models come with Popoff valve mask, oxygen tank.

Code	Description	Selling unit
39011000	RESUSCITATOR PEDIATRIC WITH SILICONE 500 ML > 6-30 KG MASK NR 2 AND 3	CX. 1
39011035	RESUSCITATOR PEDIATRIC WITH SILICONE 550 ML > 6-30 KG MASK NR 2 AND 3 AND 0,6 L BALOON	CX. 1
39011016	RESUSCITATOR NEO-NATAL WITH SILICONE 280 ML < 6 KG MASK NR 0 AND 1	CX. 1
39011039	RESUSCITATOR ADULT DISPOSABLE 1700 ML > 30 KG MASK NR 4 AND 5 AND 2,5L BALOON	CX. 1
39011004	RESUSCITATOR ADULT WITH SILICONE 1700 ML > 30 KG MASK NR 4 AND 5	CX. 1
39011034	RESUSCITATOR ADULT WITH SILICONE 1700 ML > 30 KG MASK NR 4 AND 5 AND 2,5L BALOON	CX. 1
39011001	RESUSCITATOR ADULT WITH SILICONE	CX. 1

FACTORMED
Produtos Médicos & Hospitalares