

# ENDOTRACHEAL TUBE WITH CUFF



## MAIN SPECIFICATIONS

- Atraumatic hole in the tip of the tube.
- For oral or nasal intubation.
- For a better intubation the tube has inches in depth markings which indicate the distance from the distal end.
- With radiopaque line.
- Indicated for the management of ventilation through the nasal or oral intubation of the trachea.
- With cuff.

Code	Description	Selling unit
09011058	ENDOTRACHEAL TUBE WITH CUFF NR 3,5	CX. 10
09011053	ENDOTRACHEAL TUBE WITH CUFF NR 6,5	CX. 10
09011057	ENDOTRACHEAL TUBE WITH CUFF NR 4	CX. 10
09011067	ENDOTRACHEAL TUBE WITH CUFF NR 5	CX. 10
09011056	ENDOTRACHEAL TUBE WITH CUFF NR 4,5	CX. 10
09011068	ENDOTRACHEAL TUBE WITH CUFF NR 5,5	CX. 10
09011052	ENDOTRACHEAL TUBE WITH CUFF NR 6	CX. 10

09011054	ENDOTRACHEAL TUBE WITH CUFF NR 7,5	CX. 10
09011024	ENDOTRACHEAL TUBE WITH CUFF NR 7	CX. 10
09011025	ENDOTRACHEAL TUBE WITH CUFF NR 8	CX. 10
09011046	ENDOTRACHEAL TUBE WITH CUFF NR 8,5	CX. 10
09158000	ENDOTRACHEAL TUBE WITH CUFF NR 2,5	CX. 10
09158001	ENDOTRACHEAL TUBE WITH CUFF NR 3	CX. 10
09011069	ENDOTRACHEAL TUBE WITH CUFF NR 9,5	CX. 10
09011043	ENDOTRACHEAL TUBE WITH CUFF NR 9	CX. 10
09011070	ENDOTRACHEAL TUBE WITH CUFF NR 10	CX. 10
09138001	ENDOTRACHEAL TUBE WITH CUFF NR 10,5	CX. 10
09138002	ENDOTRACHEAL TUBE WITH CUFF NR 11	CX. 10