

ENDOTRACHEAL TUBE WITHOUT CUFF



MAIN SPECIFICATIONS

- Atraumatic hole in the tip of the tube.
- For oral or nasal intubation.
- For a better intubation the tube has inches in depth markings which indicate the distance from the distal end.
- With radiopaque line.
- Indicated for the management of ventilation through the nasal or oral intubation of the trachea.
- Without cuff.

Code	Description	Selling unit
09011055	ENDOTRACHEAL TUBE WITHOUT CUFF NR 3,5	CX. 10
09011060	ENDOTRACHEAL TUBE WITHOUT CUFF NR 4	CX. 10
09011086	ENDOTRACHEAL TUBE WITHOUT CUFF NR 3	CX. 10
09011061	ENDOTRACHEAL TUBE WITHOUT CUFF NR 7	CX. 10
09011062	ENDOTRACHEAL TUBE WITHOUT CUFF NR 7,5	CX. 10