

ENDOTRACHEAL TUBE REINFORCED WITH CUFF



MAIN SPECIFICATIONS

- Endotracheal tube transparent wireframe, radiopaque, graduated bold line to mark the insertion site
- 15 mm removable connector to pass through the laryngeal mask.
- End in silicone and stabilizer bar to prevent removal of the tube through the mask.
- Cuffed.

Code	Description	Selling unit
09158002	ENDOTRACHEAL TUBE REINFORCED WITH CUFF NR 7	CX. 10
09158003	ENDOTRACHEAL TUBE REINFORCED WITH CUFF NR 7,5	CX. 10
09158004	ENDOTRACHEAL TUBE REINFORCED WITH CUFF NR 8	CX. 10