

CATHETER IV WITH FIXATION WINGS



MAIN SPECIFICATIONS

- Conventional intravenous catheter.
- With wings.
- Cannula made of polyurethane, it is latex and PVC free.
- Hydrophobic filter.
- Transparent viewing chamber.
- Single use.
- Sterilized in ethylene oxide.

Code	Description	Selling unit
09001011	CATHETER IV WITH FIXATION WINGS 14G/2'' (2,2X50 MM)	CX. 50
09001021	CATHETER IV WITH FIXATION WINGS 18GX1/4 V.DORSAL	CX. 50
09001022	CATHETER IV WITH FIXATION WINGS 20G/1 1/4'' (1,1x 32 MM)	CX. 50
09001023	CATHETER IV WITH FIXATION WINGS 22G/1'' (0,9x25 MM)	CX. 50
09001024	CATHETER IV WITH FIXATION WINGS 24G/3/4'' (0,7x19 MM)	CX. 50
09001027	CATHETER IV WITH FIXATION WINGS 18GX3/4	CX. 50
09001029	CATHETER IV WITH FIXATION WINGS 16G/2'' (1,7X50 MM)	CX. 50

FACTORMED
Produtos Médicos & Hospitalares