

# CATHETER IV WITHOUT FIXATION WINGS



## MAIN SPECIFICATIONS

- Conventional intravenous catheter.
- Without wings.
- Cannula made of polyurethane, it is latex and PVC free.
- Hydrophobic filter.
- Transparent viewing chamber.
- Single use.
- Sterilized in ethylene oxide.

Code	Description	Selling unit
09001014	CATHETER IV WITHOUT FIXATION WINGS 14G/2 (2,2X50 MM)	CX.100
09001015	CATHETER IV WITHOUT FIXATION WINGS 16G/2" (1,7X50 MM)	CX.100
09001016	CATHETER IV WITHOUT FIXATION WINGS 18G/2" (1,3X45 MM)	CX.100
09001017	CATHETER IV WITHOUT FIXATION WINGS 20G/1 1/4" (1,1X32 MM)	CX.100
09001018	CATHETER IV WITHOUT FIXATION WINGS 22G/1" (0,9X25 MM)	CX.100
09001020	CATHETER IV WITHOUT FIXATION WINGS 24GX3/4" (0,7X19 MM)	CX.100

09001025	CATHETER IV WITHOUT FIXATION WINGS 20GX1 1/4" (1,1X32 MM)	CX.100
09001049	CATHETER IV WITHOUT FIXATION WINGS 18GX1 3/4" (1,3X45 MM)	CX.100
09001038	CATHETER IV WITHOUT FIXATION WINGS 22GX1" (0,9X25 MM)	CX.100