

# SILICONE REUSABLE MASK (AMBU)



## MAIN SPECIFICATIONS

- Silicone material, latex-free.
- Easy disassembly / assembly for effective cleaning.
- Transparent mask with multifunctional coverage.
- Connection for ventilation.
- Individual packaging with instructions for use.
- Designed for a greater degree of grip during use.
- The anatomical design of the masks allows for better comfort in the face of the patient.
- Fully detachable, thus facilitating the storage and cleaning of all components.

Code	Description	Selling unit
39011027	SILICONE REUSABLE MASK NR 0 (NEONATAL)	CX. 1
39011028	SILICONE REUSABLE MASK NR 0 (PEDIATRIC S)	CX. 1
39011020	SILICONE REUSABLE MASK NR 3 (PEDIATRIC L)	CX. 1
39011019	SILICONE REUSABLE MASK NR 2 (PEDIATRIC M)	CX. 1
39011018	SILICONE REUSABLE MASK NR 4 (ADULT M)	CX. 1