

LARYNGEAL MASK AIRWAY IN SILICONE



MAIN SPECIFICATIONS

- A disposable silicone laryngeal mask airway.
- The cuff adapts to the contour of the oropharyngeal area to provide a secure seal.
- Medical grade silicone.
- Crush Resistant - The main tube wall thickness reduces the risk of crushing.
- Color Coded Pilot Balloon - Ease of identification of sizes by color coded pilot balloons.
- Secure Seal for Airway Management - Highly elastic silicone material allows the cuff of the mask to form itself to the contour of the oropharyngeal area and provide a seal.

Code	Description	Selling unit
39158010	LARYNGEAL MASK AIRWAY IN SILICONE DISPOSABLE NR 1	CX. 1
39158011	LARYNGEAL MASK AIRWAY IN SILICONE DISPOSABLE NR 2	CX. 1
39158012	LARYNGEAL MASK AIRWAY IN SILICONE DISPOSABLE NR 3	CX. 1
39158013	LARYNGEAL MASK AIRWAY IN SILICONE DISPOSABLE NR 4	CX. 1