

LARYNGEAL MASK AIRWAY I-GEL



MAIN SPECIFICATIONS

- Single use.
- Establishing an anatomical seal, not inflatable pharyngeal structures and laryngeal perilaryngeal, while avoiding compression trauma.
- With thermoplastic elastomer medical grade, which in contact with body heat, swells making the seal between the mask and the hypopharynx.
- Available in three sizes for adult and with color code, each mask is involved in a colored code case.

Code	Description	Selling unit
39100005	LARYNGEAL MASK I-GEL NR 3 YELLOW (30-60KG)	CX. 10
39100006	LARYNGEAL MASK I-GEL NR 4 GREEN (50-90KG)	CX. 10
39100007	LARYNGEAL MASK I-GEL NR 5 ORANGE (+90KG)	CX. 10