

OXYGEN TUBE



MAIN SPECIFICATIONS

- Transparent, flexible tube (without forming folds or creases that prevent the passage way of medical gases) with two connectors that allow you to adjust to most adapters.
- Latex free.
- Used for connection to medical gas cylinders or adapters.
- Available in different sizes.

Code	Description	Selling unit
39011031	OXYGEN TUBE 210 CM WITH 2 CONNECTORS	CX. 1
39011032	OXYGEN TUBE (3MM OF DIAMETER) 30 M ROLL (GREEN)	CX. 1
39016027	OXYGEN TUBE 15 M WITH 2 CONNECTORS	CX. 1
39037005	OXYGEN TUBE (7MM DIAMETER) 50 M ROLL	CX. 1