

CERVICAL COLLAR



MAIN SPECIFICATIONS

- Moldable. Comfortable. Anti-allergy. Washable.
- Semi-rigid cervical immobilization.
- Soft foam with coverage in 88% cotton and 12% polyurethane.
- Anti-allergy. Washable. Adjustable with velcro closure.
- Its anatomical shape allows a perfect fit in the spine.
- For myalgias the cervical area.
- Hard "Schanz".
- Rigid cervical immobilization.
- Made from polyethylene vented to air circulation.
- Padded shoulder to chin.
- Adjustable with velcro closure.
- The anatomical shape allows a better therapeutic support with less discomfort to the user.
- Secure fastening system.
- For sprains trauma and myalgias cervical area.

Code	Description	Selling unit
39080003	CERVICAL COLLAR XS (GREEN)	CX. 1
39080004	CERVICAL COLLAR M (ORANGE)	CX. 1
39080005	CERVICAL COLLAR L (BLUE)	CX. 1
39080012	CERVICAL COLLAR XL (DARKER BLUE) STOUT	CX. 1
39080016	CERVICAL COLLAR (HEIGHT ADJUSTABLE 4 POSITIONS)	CX. 1
39080000	CERVICAL COLLAR S (YELLOW)	CX. 1