

## AEROSOL MASK WITH CUP



### MAIN SPECIFICATIONS

- Aerosol mask with a nebulizer cup and 6 ml capacity tube for oxygen combined therapy with drug and humidification.
- Single use.
- Latex free.
- Available for child and adult.

Code	Description	Selling unit
39279002	ADULT AEROSOL MASK WITH NEBULIZER CUP 6 ML	CX.100
39279003	PEDIATRIC AEROSOL MASK WITH NEBULIZER CUP 6 ML	CX.100