

POCKET MASK AND FILTER FOR THE MASK



MAIN SPECIFICATIONS

- Suitable for adults, children and infants.
- The mask is reusable and easy disinfection.
- Transparent material for better visualization of the procedure.
- Removable check valve of high efficiency.
- Soft cushion to prevent air leakage and enable greater patient comfort.
- Valve disposable to prevent cross contamination.
- Elastic to keep the mask attached to the patient.
- 15 mm connector for coupling of respiratory support equipment, such as resuscitators and oxygen.
- Latex-free material and nontoxic.
- It is supplied ready for first use.
- Supplied in a compact plastic case for mask storage, with a corrugated tube for an easier use.

Code	Description	Selling unit
39011002	POCKET MASK	CX. 1

FACTORMED
Produtos Médicos & Hospitalares