

DUMMY NEWBORN FOR BASIC LIFE SUPPORT



MAIN SPECIFICATIONS

- Faithful reproductions of natural anatomical features and weight, which can be achieved by filling the manikin with water.
- The simulation is realistic chest expansion, as well as cardiac compression achieved with a semi-rigid thorax.
- Thanks to an advanced electrical pad for visible warnings, shows that the resuscitation procedure is correct regarding the ventilation of the lungs, super compression of the chest, exact position of the fingers on the sternum and cardiac massage.
- Is provided with a suitable plastic traveling bag, which can also be used as a cushion during training.

Code	Description	Selling unit
39080023	DUMMY NEWBORN FOR BLS	CX. 1