

## STRECHER BELTS KIT



### MAIN SPECIFICATIONS

- Light quick release buckle, intuitive and safe long term.
- The major components of the buckle ensure mechanical and chemical characteristics
- Great technical strength.
- Stiffness optimized.
- Enhanced dimensional stability in hot conditions.
- Maximum resistance to temperatures well below freezing.
- Good fire resistance.
- Good resistance to oil, grease and solvents.
- Width: 65 mm
- Weight: 50 g
- Load weight limit: 50 kg
- Set of three belts with carrying case.

Code	Description	Selling unit
39080028	STRECHER BELTS KIT	CX. 3