

LONG HARD PLAN



MAIN SPECIFICATIONS

- Plan long which reduces the possibility of contamination during transport of the patient.
- Fully recyclable high-density polyethylene molded in one piece.
- With rigidity and light weight characteristics.
- Resistance to deformation and corrosion.
- Fully waterproof.
- It does not absorb liquid.
- Immune to infiltration.
- Easy to clean (water and soap are sufficient).
- Fully compliant radio (the patient can be removed even after the litter ray or CT).

Code	Description	Selling unit
39080006	LONG HARD PLAN	CX. 1