

## HEAD IMMOBILIZER



### MAIN SPECIFICATIONS

- The head immobilizer ensures optimum immobilization of head injured patients.
- When used in conjunction with a column immobilized on a stretcher and with the hard cervical collar, allows the transport of a patient in critical condition during long trips or on uneven.
- Provides support to the skull, avoiding as much as possible a greater compression of the skull and completing the immobilization of the spine through the cervical collar.
- Consisting of two mono blocks, made of soft plastic material and having a base with integrated belts adaptable to various types of stretchers.
- The two shells are positioned in a special base through a very strong fastening system in Velcro which may be fixed in any type of rigid litter supplied through the system or using Velcro belts. Both shells have deliberately made holes to permit inspection of the external ear, can verify the occurrence of any loss of blood or fluid, in addition to allowing rescue operators communicate with the patient.
- The database allows the setting of only two blocks for adult or pediatric patients, due to the additional pad can be removed.
- The material used in the production of the blocks is waterproof, long term, preventing the absorption of body fluids (blood, vomit, mucus).
- It can be damaged by chemicals or solvents, or altered by physical agents and remains soft even when exposed to temperature variations.

- It is easily washable and is radio-transparent.

| Code     | Description      | Selling unit |
|----------|------------------|--------------|
| 39080007 | HEAD IMMOBILIZER | CX. 1        |