

ASKINA TRANSORBENT



MAIN SPECIFICATIONS

- The Transorbent dry hydrogel technology features a patented design that provides a unique absorption process. It is able to transfer fluid away from the wound, capture it, and eliminate the excess moisture.
- Superior absorption capacity.
- Does not stick to the wound surface.
- Cushioning effect due to the polyurethane foam.
- No maceration or residue.
- Can be cut to individually adapted sizes and shapes.
- Pressure ulcers, venous leg ulcers, arterial leg ulcers, diabetic foot ulcers, traumatic wounds, surgical wounds, 1st and 2nd degree burns, aid in the prevention of skin breakdown.

Code	Description	Selling unit
01001000	ASKINA TRANSORBENT 10X10CM	CX. 5
01001005	ASKINA TRANSORBENT 5X7CM	CX. 10
0072790V	ASKINA TRANSORBENT 15X15CM	CX. 5
0072791W	ASKINA TRANSORBENT 20X20CM	CX. 5
0072710C	ASKINA TRANSORBENT BORDER 14X14CM	CX. 5

0072715D	ASKINA TRANSORBENT BORDER 17X17CM	CX. 5
0072720E	ASKINA TRANSORBENT BORDER 24X24 CM	CX. 5
0072705B	ASKINA TRANSORBENT BORDER 9X14CM	CX. 5
0072792X	ASKINA TRANSORBENT SACRUM 18X20CM	CX. 5