

ASKINA GEL HYDROGEL



MAIN SPECIFICATIONS

- Askina Gel is a clear, viscous, sterile hydrogel. Dependant on wound conditions, Askina Gel rehydrates dry necrotic tissue, to promote debridement in dry wounds, or absorbs wound slough and exudate where a small amount of exudate is present.
- Askina Gel maintains a moist wound healing environment conducive to natural healing.
- Easy to apply, in any position, with one hand.
- Remains in place after application.
- A special cannula allows to reach deeper wounds sinuses, fistulas.
- More hydration: donates 40% of its water to the wound
- Crystal clear, provides a more accurate wound assessment
- Low exuding, sloughy or dry necrotic wounds: venous and arterial leg ulcers, pressure ulcers, diabetic foot ulcers and traumatic wounds.
- Contraindication: Do not use on patients with known sensitivity to the gel or one of its components; Do not use on third degree burns; Do not use if tube is damaged or seal is broken.

Code	Description	Selling unit
01001001	ASKINA GEL 15 GRS	CX. 5

FACTORMED
Produtos Médicos & Hospitalares