

# ASKINA BIOFILM TRANSPARENT



## MAIN SPECIFICATIONS

- Askina Biofilm transparent is a thin sterile transparent hydrocolloid dressing for moist wound treatment with the capacity to absorb residual exudate.
- Composition: Sodium-carboxymethylcellulose, Polyisobutylene and Polyurethane film.
- Thin design that fits any body contours.
- Transparent, allows visual control of the wound.
- Easy to apply.
- Stays in place during showering.
- Management of low exuding superficial chronic wounds (pressure ulcers, leg ulcers), superficial acute wounds (1<sup>st</sup> and 2<sup>nd</sup> degree burns, abrasions) and post operative wounds.

Code	Description	Selling unit
01001003	ASKINA BIOFILM TRANSPARENT 10X10CM	CX. 10
F72092	ASKINA BIOFILM TRANSPARENT 20X20CM	CX. 5
F72091	ASKINA BIOFILM TRANSPARENT 15X15CM	CX. 5
F72095	ASKINA BIOFILM TRANSPARENT 5X20 CM	CX. 10