

ASKINA CARBOSORB CHARCOAL DRESSING



MAIN SPECIFICATIONS

- Askina Carbosorb is a conformable, sterile wound dressing with a thick, absorbent layer of non woven viscose-rayon, and an activated charcoal cloth layer for bacteria and malodorous absorption.
- Askina Carbosorb has limited absorption capacity: it can be used alone on low exuding malodorous wound, or in combination with other absorbing dressings for moderately to heavily exuding malodorous wounds.
- Efficient odour absorption.
- Can be cut.
- Askina Carbosorb may be used for the management of low to heavily exuding, partial to full thickness chronic wounds, in particular when infection and/or offensive odour is present: venous leg ulcers, arterial leg ulcers, pressure ulcers, fungating wounds, infected and/or malodorous surgical or traumatological wounds.

Code	Description	Selling unit
01001004	ASKINA CARBOSORB CHARCOAL DRESSING 10X10CM	CX. 10
01001015	ASKINA CARBOSORB CHARCOAL DRESSING 10X20CM	CX. 10

FACTORMED
Produtos Médicos & Hospitalares