

ACTISORB SILVER SILVER AND CHARCOAL DRESSING



MAIN SPECIFICATIONS

- Actisorb Silver 220 Activated Charcoal Dressing with Silver is an activated charcoal dressing with silver, enclosed in a non-adherent nylon sleeve.
- Actisorb Silver 220 is different from other charcoal dressings, it has a triple action that not only masks odour but also addresses its cause.
- Antimicrobial action - Actisorb Silver 220 is effective against 150 microorganisms associated with wound infection, including MRSA and VRE.
- A recent study shows the dressing provides effective antimicrobial action without causing any harm to host cells nor preventing fibroblast proliferation.
- Odour and exudate management - Actisorb Silver 220 eliminates malodour thereby improving the patient's quality of life. As malodour might be accompanied by wound fluids, Actisorb Silver 220 is able to manage exudate and helps maintain a moist wound healing environment.
- Actisorb Silver 220 is indicated to manage bacterial contamination, infection or odour in a wide range of wounds, including: Pressure ulcers; Leg ulcers deriving from different aetiologies; Diabetic foot ulcers; Fungating carcinoma; Traumatic and surgical wounds where bacterial contamination, infection or odour occurs.

Code	Description	Selling unit
01108004	ACTISORB SILVER 220 6,5X9,5 CM	CX. 10
01108005	ACTISORB SILVER 220 10,5X10,5 CM	CX. 10
01108006	ACTISORB SILVER 220 19X10,5 CM	CX. 10