

# OMNIFILM MICRO POLYACRYLATE ADHESIVE



## MAIN SPECIFICATIONS

- Hypoallergenic adhesive on porous plastic film with adhesive polyacrylate. Permeable to air and steam, impervious to dirt, water resistant and with good adhesion.
- It removed painlessly and without leaving marks.
- Easy to cut.
- Radiolucent.
- Adhesive coverage and retention.
- Transparent fixation for probes, cannulas, etc.
- For sensitive skin.

Code	Description	Selling unit
03046000	OMNIFILM ADHESIVE 5 M X2,5 CM	CX. 1
03046001	OMNIFILM ADHESIVE 9,2 M X1,25 CM	CX.1
03046002	OMNIFILM ADHESIVE 9,2 M X 2,5 CM	CX. 1
03046022	OMNIFILM ADHESIVE 5 M X 1,25 CM	CX. 1

03046023	OMNIFILM ADHESIVE 5 M X 5 CM	CX. 1
----------	------------------------------	-------

---

03046025	OMNIFILM ADHESIVE 9,2 M X 5 CM	CX. 1
----------	--------------------------------	-------

---