

# PARANET PARAFFIN GAUZE STERILE



## MAIN SPECIFICATIONS

- Paranet sterile gauze made of 100% cotton open loop Leno weave impregnated with paraffin, yellow soft (has a supersaturation of the solution) for use in superficial wounds, burns, skin grafts and traumatic injuries.
- It is recommended as a contact layer with the primary wound.
- Paranet be used in conjunction with an appropriate secondary dressing.
- Paranet is sterile gauze made of 100%. Presented in gauze and bandages multi-piece strips Paranet is packaged in polypropylene containers and sealed with a tamper proof seal in individual aluminum bag.

Code	Description	Selling unit
01016009	PARANET PARAFFIN GAUZE STERILE 10X10 CM PACKED INDIVIDUALLY	CX. 10
01016050	PARANET PARAFFIN GAUZE STERILE 10CMX7M	CX. 1