

# TIELLE PLUS HYDROPOLYMER DRESSING



## MAIN SPECIFICATIONS

- Tielle Plus Hydropolymer Adhesive Dressing is an adhesive hydropolymer foam dressing with LiquaLock Advanced Absorption Technology for medium to highly exuding wounds.
- Tielle Plus helps maintain a moist wound healing environment through the combination of absorbency and moisture vapour transfer.
- The hydropolymer with LiquaLock Technology: Expands and conforms to the wound bed, ensures fluid is locked away, minimising the risk of leakage and maceration.
- Tielle Plus dressings contribute to a decrease in: Maceration, Excoriation, Leakage, Number of dressing changes, Skin reactions (to the adhesive border).
- The unique adhesive allows: Repositioning as needed on initial application and dressing flexibility to conform to difficult to dress areas.
- Tielle Plus can be used in medium to highly exuding, partial and full thickness wounds, including: Pressure ulcers; Lower extremity ulcers, (e.g. venous, arterial and mixed aetiology), Diabetic foot ulcers and Donor sites.

Code	Description	Selling unit
01108000	TIELLE PLUS HYDROPOLYMER DRESSING 11X11 CM	CX. 10