

AQUACEL EXTRA HYDROFIBER DRESSING



MAIN SPECIFICATIONS

- Aquacel Extra wound dressing is composed of two layers of Hydrofiber Technology stitched together. It is nine times stronger and has 39% greater absorbency as compared to the original Aquacel wound dressing.
- Aquacel Extra Hydrofiber dressing is specifically suited to help manage moderate to highly exuding wounds and is cleared for the same indications as the original Aquacel Extra dressing. Aquacel Extra wound dressing is constructed with two-dimensional strengthening fibers that are stitched into the dressing. While made from the same base material that is used to manufacture Hydrofiber Technology, the fibers do not gel on contact with fluid so that they maintain their integrity and become encompassed by the gel, providing additional strength and aiding the easy removal.
- Dressings with Hydrofiber Technology gel on contact with fluid, providing clinicians and patients with substantial wound care benefits. This unique gelling action enables the dressings to lock in exudate and its harmful components and contour closely to the wound bed, in response to changing wound conditions. Aquacel dressings are soft, absorbent non-woven wound dressings based on Hydrofiber Technology.

Code	Description	Selling unit
01101017	AQUACEL EXTRA DRESSING 5X5 CM	CX. 10
01101013	AQUACEL EXTRA DRESSING 10X10 CM	CX. 10
01101018	AQUACEL EXTRA DRESSING 15X15 CM	CX. 5

FACTORMED
Produtos Médicos & Hospitalares