

# VARIHESIVE EXTRA MINCE HYDROCOLLOID DRESSING



## MAIN SPECIFICATIONS

- Vahiresive Extra Mince is a hydrocolloid dressing for local treatment of pressure ulcers and slightly exudative lesions.
- Consists of an inner layer of hydrocolloid (in contact with the lesion) contained within a polymeric adhesive matrix.
- The dressing provides a moist environment for the wound which promotes healing, aiding the autolytic debridement, facilitating the migration of epithelial cells, and allowing the removal of non-traumatic dressing without damaging the newly formed tissue.
- The adhesive layer containing elastomeric polymers that increase the ability of the dressing to handle the wound drainage by forming a cohesive gel.
- Vahiresive Extra Mince is a dressing that its composition and design, is particularly suitable for the treatment of exudative skin lesions slightly. It is indicated for pressure ulcers and vascular ulcers.
- Its use is contraindicated in patients with known sensitivity dressing or its components
- Available in different sizes.

Code

Description

Selling unit

01101004	VARIHESIVE EXTRA MINCE 7,5X7,5 CM	CX. 10
01101005	VARIHESIVE EXTRA MINCE 10X10 CM	CX. 10
01101006	VARIHESIVE EXTRA MINCE 15X15 CM	CX. 10