

## OMNISTRIP CLOSURE STRIP



### MAIN SPECIFICATIONS

- Omnistrip are wound closure strips made of skin coloured non-woven fabric.
- The low-irritant properties of Omnistrip are enhanced by the high permeability to air and water vapour of the material.
- Flexible, non-woven material allows for expansion arising from oedema and flexion.
- Acrylic adhesive reduces the possibility of allergic reaction.
- Rounded edges help prevent roll up.
- Does not absorb X-rays and may be left in place for radiological exposure.
- Cutaneous suture strips that ensure a reliable and direct closing small incisions and primary as well as the alleviation of wound sutures.
- The non-woven with hypoallergenic adhesive skin color has an excellent tensile strength.
- Allows a secure fixing of the wound edges and can be removed without causing pain.
- Rounded corners for a good grip.
- Permeable to air and water vapor.
- Sterile.

Code	Description	Selling unit
------	-------------	--------------

01046010	OMNISTRIP 6X76 MM	CX. 150
01046011	OMNISTRIP 3X76 MM	CX. 250
01046012	OMNISTRIP 6X38 MM	CX. 300
01046013	OMNISTRIP 12X101 MM	CX. 300
01046043	OMNISTRIP 6X101 MM	CX. 500