

# LEUKOSAN STRIP CLOSURE STRIP



## MAIN SPECIFICATIONS

- Leukosan strip are wound closure strips made of skin coloured non-woven fabric.
- The low-irritant properties of Leukosan strip are enhanced by the high permeability to air and water vapour of the material.
- Flexible, non-woven material allows for expansion arising from oedema and flexion.
- Acrylic adhesive reduces the possibility of allergic reaction.
- Rounded edges help prevent roll up.
- Does not absorb X-rays and may be left in place for radiological exposure.

Code	Description	Selling unit
01254021	LEUKOSAN STRIP 6X75 MM	CX.150
01254020	LEUKOSAN STRIP 3X76 MM	CX. 250
01254014	LEUKOSAN STRIP 6X38 MM	CX. 300
01254023	LEUKOSAN STRIP 12X101 MM	CX. 300
01254022	LEUKOSAN STRIP 6X100 MM	CX. 500