

HYDROCOLL THIN HYDROCOLOID DRESSING



MAIN SPECIFICATIONS

- Hydrocolloid dressing is a sterile, self adhesive hydrocolloid wound dressing that features a bevelled edge for secure adhesion to surrounding skin.
- Hydrocoll forms a soft gel when in contact with wound fluid, producing a moist wound environment suitable for autolytic debridement and protection of granulating tissue.
- Hydrocoll at a glance: Absorbs and retains exudate, while maintaining an adequate moisture balance.
- Thin, flexible material for conformability. Bacteria and water proof.
- Bevelled edge for secure adhesion.
- Low profile resists roll-off.
- Can be removed painlessly and without leaving residue.
- Available in specifically designed shapes for problem-free application in difficult areas.
- Can remain on a wound for up to 5 days.

Code	Description	Selling unit
01046031	HYDROCOLL THIN DRESSING 7,5X7,5 CM	CX. 10

